

Bonaduz, July 27, 2017

Reprocessing of Hamilton Medical consumables

High-level disinfection of Hamilton Medical flow sensor and expiratory valve sets

Compatibility with STERIS high-level disinfectants

The Hamilton Medical autoclavable expiratory valve sets and autoclavable flow sensor are suitable for high-level disinfection with the following STERIS high-level disinfectants:

- Revital-Ox Resert High-Level Disinfectant
- Revital-Ox Resert High-Level Disinfectant – Chemosterilant (5-minute cycles)
- Resert XL HLD High-Level Disinfectant (8-minute cycles)

The following consumables may be successfully high-level disinfected for up to 50 cycles:

PN161188: Expiratory valve set for HAMILTON-C1/T1/MR1, neonatal, autoclavable
PN161175: Expiratory valve set for HAMILTON-C1/T1/MR1, adult/pediatric, autoclavable
PN151228: Expiratory valve housing for HAMILTON-G5/S1, all patient groups, autoclavable
PN151233: Expiratory valve membrane for HAMILTON-G5/S1, all patient groups, autoclavable
PN160245: Expiratory valve membrane for HAMILTON-C2/C3/C6, all patient groups, autoclavable
PN950185: Flow sensor, adult/pediatric, autoclavable

With kind regards

Hamilton Medical AG



Stefanie Schmidt

Product- and Hygienemanager
stschmidt@hamilton-medical.com

☎ +41 58 610 20 88

📠 +41 58 610 00 20



Ray Curtis

Quality Engineer
rcurtis@hamilton-medical.com

☎ +41 58 610 21 37

📠 +41 58 610 00 20