


Stretchers

Reliability, safety and ease: a better outcome



Enhancing outcomes for patients and their caregivers:

Hill-Rom

The Hill-Rom[®] mobility is life[™] concept

Hill-Rom is committed to the concept that mobilising patients early and supporting their independence improves life – in the hospital and at home.

- By enabling caregivers to provide efficient care, have more time at the bedside and get patients mobilised safely, positive outcomes can be achieved. Easier and more pleasant working conditions for caregivers allow more time for the dignified care of patients, which can **accelerate recovery**.
- Achieving this requires a **commitment to continuous innovation**, assistive technology, quality and clinical validation.
- Hill-Rom thinks, develops and manufactures based on this and always with the mission of **enhancing outcomes for patients and caregivers**.



Stretchers, the **unsung heroes** of the healthcare environment

They take a beating in the toughest conditions, but still do their job reliably and without fail. Without seeking the spotlight, stretchers play a pivotal role in helping you deliver the best patient care.

At Hill-Rom, we believe in the importance of stretchers. But we see reliability as only part of the equation, making a difference in caregiver safety and ease is just as critical. We've designed our family of stretchers to meet unique caregiver needs across a variety of settings, from procedural and transport to OB-GYN and surgical.

The result, a safe, easy stretcher experience - a better outcome for everyone. With stretchers from Hill-Rom, you'll see the results and feel the difference.

Feel the difference

A stretcher from Hill-Rom is more than a smart investment. It's a deeper connection with a dedicated long-term partner.

Every Hill-Rom stretcher comes with the reliability and best-in-class service you expect from us. We're committed to helping achieve a better outcome for patients and providers at every point of care, and our family of stretchers is another great example.

Capacity

Hill-Rom is the only hospital stretcher company to offer a 318 kg (700 lbs) weight capacity for all of its stretchers. This versatility makes Hill-Rom stretchers perfect emergency stretchers, able to handle almost any circumstance.

Safety

Our **BackSaver Fowler™** Feature uses the patient's own weight to help you raise and lower the head section of the stretcher easily and comfortably.

Service

Our hospital stretcher service warranty is supported by service centres and Hill-Rom certified engineers located throughout the country.

Stability

An industry first, our **Auto Contour™** Feature allows automatic rising of the knee section as the fowler is raised. This helps minimise patients slipping down the mattress and minimises patient repositioning.

Control

A convenient hand brake helps you negotiate ramps and corners without using your body weight, resulting in safe and easy stops.

Accuracy

Our highly accurate built-in scale displays both kilograms and pounds simultaneously, eliminating the need to toggle between the two. The built-in scale is only available in countries which do not require a CE mark.



The industry's most extensive stretcher portfolio



Surgical Stretcher*

- Flexible head section
- Surgical rail for accessories



Transport Stretcher

- Optional fifth wheel mobility
- Ergonomic push handles

*Whatever your clinical
need, we have a stretcher
to meet it and help you
achieve a better outcome
in the process.*



Trauma Stretcher

- Full length X-Ray cassette holder
- X-Ray lift handles on both sides

4-Sided Brake and Steer
Optional on all stretchers



Electric Stretcher*

- Patient head and knee controls
- Lockout controls for caregivers



Procedural Stretcher

- Tuckaway siderails for zero transfer gap
- Fifth wheel mobility
- BackSaver Fowler
- Auto Contour



OB-GYN Stretcher*

- Integrated foot supports
- Stowable foot section

*The difference is in
the details.*

Every stretcher in our portfolio contains special features to ensure a better outcome whenever possible. The industry leading capacity of all Hill-Rom stretchers is 318 kg (700 lbs).

Hill-Rom is a leading global medical technology company with more than 6,500 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

mobility is life 

Hill-Rom is committed to the concept that mobilising people early and supporting their independence, improves life - in the hospital and at home.

This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.



These medical devices bear the CE mark and comply with all applicable regulations. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Medical beds manufacturer:

Hill-Rom SAS - B.P. 14 - Z.I du Talhouët - 56330 Pluvigner - France (also European Representative)

Völker GmbH - Wullener Feld 79 - 58454 Witten - Germany

Hill-Rom Holdings, Inc. - 1069 State Route 46 East - Batesville, IN 47006 - USA

Lifts manufacturer: Liko AB - Nedre vägen 100 - 975 92 Luleå - Sweden

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2014 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN145308-01, 18-07-2014

Not all products/options are available in all countries.

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

Enhancing outcomes for
patients and their caregivers:

Hill-Rom